BODY FLEX®

GUN 5

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ENGLISH

Attachment heads

- U-shaped head - for the muscles alongside the spine and around the neck; do not apply pressure directly to the vertebrae.
- Large round head - for large muscle groups such as thighs, glutes, calves and back.
- Small round head - for joints, trigger points and smaller treatment areas.
- Small flat head - for general use on most parts of the body.

Operation & charging

- Press and hold the Start/Stop button to turn the device ON/OFF.
- Press the button briefly to change the speed level.
- Pull the attachment straight out to remove it; push the new one in firmly.
- Charge the device with a 5V USB power source. Before first use, charge for about 3 hours.
- The indicator lights show the charging status, battery level and selected speed.

Safe use

- Use on clean, dry skin with light to moderate pressure.
- Treat one area for about 15-60 seconds, then move to the next area.
- Use only on soft tissue. Avoid bones, the head, open wounds and swollen or inflamed areas.
- Mild redness may occur; stop immediately if you feel sharp pain, numbness or discomfort.
- Keep fingers, hair and loose clothing away from the moving parts and air vents.
- Keep the device dry. Do not immerse in water.

Do not use / consult a doctor first

- If you have a pacemaker, defibrillator or serious heart condition.
- If you have a bleeding disorder, take anticoagulants, or have recently had surgery or a fracture.
- If you have acute inflammation, tumors, are pregnant, or are undergoing specialist rehabilitation.
- Do not use on the head, face, cervical spine, vertebrae, or areas with very little muscle cover.
- Keep the device out of reach of children and never leave it unattended while charging.

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